**METTLER TOLEDO**

## **InPro® 8000 Series**

**Turbidity sensors  
Trübungssensoren  
Sondes de turbidité**

**Instruction manual  
Bedienungsanleitung  
Instructions d'utilisation**

**English**

**page 3**

**German**

**page 21**

**French**

**page 39**

# **InPro® 8000 Series Turbidity sensors**

## **Instruction manual**

Subject to technical changes without prior notice.

© It is strictly forbidden to reprint this instruction manual or any parts thereof without the written permission of Mettler-Toledo GmbH, Process Analytics, Industrie Nord, 8902 Urdorf, Switzerland. No section or excerpt whatsoever may be reproduced or with the assistance of electronic systems be edited, duplicated or distributed, in particular in the form of photocopies, photographs, magnetic media or other recording methods. All rights reserved, especially the right of duplication and translation as well as in regard to patent and registration rights.

Mettler-Toledo GmbH, 8606 Greifensee, Switzerland

## Contents

<b>1</b>	<b>Product description</b>	<b>5</b>
1.1	Introduction	5
1.2	Equipment and scope of delivery	6
1.3	Technical data	6/7/8
<b>2</b>	<b>Safety</b>	<b>9</b>
2.1	Application compatibility	9
2.2	Proper utilization	9
2.3	Safety measures	9
2.4	Ex-Declaration	10
2.4.1	Use in Ex-zones	10
2.4.2	Ex-Classification	11
<b>3</b>	<b>Initial start-up</b>	<b>12</b>
3.1	Handling	12
3.2	Installation in vessels/reactors	12
3.3	Installation in pipes	13
3.4	Calibration in open laboratory vessels	13
3.5	Fiber optic cable extensions	13
<b>4</b>	<b>Maintenance and troubleshooting</b>	<b>14</b>
4.1	Conditions of warranty	14
4.2	Maintenance	14
<b>5</b>	<b>Removal from operation, storage, disposal</b>	<b>15</b>
5.1	Removal from operation	15
5.2	Storage	15
5.3	Disposal	15
<b>6</b>	<b>Dimensional drawings</b>	<b>16</b>
<b>7</b>	<b>Ordering information</b>	<b>17</b>
7.1	Sensors	17
7.2	Accessories	17/18
<b>8</b>	<b>Certificates</b>	<b>19</b>
8.1	Material certificate	19

# 1. Product description

## 1.1 Introduction

We thank you for your decision to purchase a METTLER TOLEDO InPro®8000 Series fiber optic measuring sensor with operation based on the principle of backward scattered light.

Please keep this instruction manual in a safe place where it is readily to hand for users of this sensor. If you have any questions regarding function of the sensor and ancillary equipment, or if you consider that the descriptions/instructions do not tally or are lacking, please contact your nearest METTLER TOLEDO Authorized Distributor, who will be pleased to provide assistance.

### Markings

The markings labeled on each sensor state the following information:

METTLER TOLEDO

Serial no. xxxxxx

Order no. 52 xxx xxx

The code following the model designation of the InPro®8100 and InPro®8200 Sensors as printed on the package label provides information on the shaft material, sealing material for the optical window (if present), and the length of the shaft in mm.

e.g. InPro®8200/S/Kalrez®-FDA/205

with S = stainless steel shaft DIN 1.4435 or 1.4404 (316L)

Kalrez®-FDA 6230 window seal

shaft length = 205 mm (8.07").

### In addition:

H = Hastelloy C276 shaft

Epoxy = Epoxy bonded window

See also chapter 7.1: Ordering information sensors.

## 1.2 Equipment and scope of delivery

### Scope of delivery

The sensors of the InPro® 8000 Series are delivered ready for use including fiber optic cable. Each sensor is accompanied by this instruction manual.

A METTLER TOLEDO Process System for the measurement of turbidity or of suspended particles is completed by a METTLER TOLEDO Turbidity Transmitter Type Trb 8300 (order no. 52 800 204). Detailed information for calibration of the measuring system can be found in the following relevant transmitter instruction manuals:

German	Order no. 52 800 243
English	Order no. 52 800 244
French	Order no. 52 800 245

### Packaging

The packaging consists of cardboard and foamed plastic. Safekeep the packaging for later use during storage or transportation of the sensor. Should you wish to dispose of the packaging material, please observe your local regulations as well as those data and instructions given in **chapter 5.3** of this manual.

### Unpacking and inspection

Please check the sensors immediately during unpacking in order to determine possible damage or missing items. Any irregularities should immediately be notified to your carrier and to your supplier.


## 1.3 Technical data

### Function

Sensors of the InPro® 8000 Series are optical sensors which measure the concentration of suspended particles in a liquid medium on the basis of the principle of backscattered light. The Transmitter Trb 8300 is equipped with a LED which beams a light – almost in the infra-red range (wavelength 880 nm) – via a fiber optic cable into the liquid medium. If the light hits on particles, then it is scattered in all directions. The light scattered at an angle of

180° (backscattered light) is captured and led back via a fiber optic cable to a photodiode in the transmitter, processed as photocurrent and the signal transformed into a value for the transmitter display.

The sensors are divided into two types, single and dual optical fiber sensors. The difference relates to the number of optical fibers present at the end (nose) of the sensor shaft.

At low particle concentrations, the measuring range is expressed in FTU (Formazin Turbidity Units). At high concentrations, the measurements are given in grams of suspended solids per liter of sample volume, using diatomaceous earth as reference substance.

Technical data

Technical data	InPro® 8200/S(H)	InPro® 8100/S	InPro® 8050
Meas. principle	dual fiber	single fiber	single fiber
Shaft material	DIN 1.4435 or 1.4404 (316L) or Hastelloy C276 (available only with Epoxy sealed window)	1.4435 or 1.4404 (316L)	PSU (Polysulfone)
Shaft diameter	12 mm	12 mm	12 mm
Shaft length	120, 205, 297 or 407 mm	120, 205, 297 or 407 mm	120 mm
Window material	Sapphire	without window	Sapphire spigot
O-ring sealing	Viton®-FDA	Viton®-FDA	Viton®-FDA
Process adaption	Pg 13.5	Pg 13.5	Pg 13.5
Plug	2 x SMA	2 x SMA	2 x SMA
Measuring range	5 FTU...4000 FTU 0...30 g/l	10 FTU...4000 FTU 0...250 g/l	10 FTU...4000 FTU 0...250 g/l
Pressure range	0...6 bar (0...87 psi) with Kalrez®-FDA sealed window 0...16 bar (0...232 psi) with Epoxy bonded window	0...6 bar (0...87 psi)	0...2 bar (0...29 psi)
Temperature range	-10...130 °C (14...266 °F) with Kalrez®-FDA sealed window -30...130 °C (-22...266 °F) with Epoxy bonded window	-30...130 °C (-22...266 °F)	0...60 °C (32...140 °F)
Steam-sterilizable	130 °C (only with Kalrez®-FDA sealed window)	130 °C	non-sterilizable

Subject to technical changes.

All sensors fall into the category described in Article 3, Paragraph 3 of the Pressure Equipment Directive, 97/23/EC.

## 2. Safety

### 2.1 Application compatibility

The wetted material parts of the sensor (shaft, window, window sealing material, O-ring) can under circumstance be non-compatible with the particular composition of the process medium and/or of the operating conditions. Responsibility to verify application compatibility lies wholly with the user.

The compatibility of different types of material are outlined on <http://www.coleparmer.com/techinfo/chemcomp.asp>. Mettler-Toledo GmbH, Process Analytics accepts no responsibility whatsoever for the correctness or accuracy of such details.

### 2.2 Proper utilization

METTLER TOLEDO InPro®8000 Series sensors are intended solely for the measurement of suspended particles or turbidity in aqueous solutions in industrial applications.

Any other use, or any operation over and above that intended by the manufacturer, are deemed to be non-permissible and incorrect, and can lead to harm or injury to material/equipment and persons. This is also relevant for applications which do not comply with the technical data of the sensor. For any damage possibly arising from such misuse, the user assumes full and sole responsibility.

### 2.3 Safety measures

The sensors InPro®8000 Series have been manufactured in line with state-of-the-art technology and in accordance with accepted technical safety regulations. Nevertheless, the sensors can still represent a source of risk and danger:

- if the sensors are operated by insufficiently trained personnel,
- if the sensors are employed incorrectly or not as intended by the manufacturer
- if the sensors are not regularly maintained or serviced.

Local legislation and regulations must be observed at all times. Such stipulations do not form an integral part of this instruction manual.

It is on principle necessary for persons handling or using the sensors to wear personal safety outfit such as protective goggles and protective clothing.

The user is responsible for the instruction and training of personnel. In this respect, additional copies of the instruction manual can be ordered from your supplier. This instruction manual is an essential element of the sensor equipment and must at all times be readily to hand for operators directly at the location of employment of the sensors.

Before the sensor is removed from the process/process adapter, it must be ensured that the process pressure has been reduced to a safe level and the process temperature lowered to a safe range. Any escape of hot process medium under pressure can cause damage to material/equipment or injury to persons.

No modification whatsoever may be carried out on the sensors. Any unauthorized modification or manipulation of the sensors results in immediate expiry of the full scope of warranty granted by the manufacturer.

## 2.4 Ex-Declaration

### 2.4.1 Use in Ex-zones

**Note:** Before installation in Ex-zones, please read and observe the following guidelines:

**Ex-Classification:**

Ⓔ II 1/2G T6/T5/T4

**Designation and number of the test certificate:**

SNCH 03 ATEX 3565 X

## 2.4.2 Ex-Classification / II 1/2G T6/T5/T4 according to EU type test certificate SNCH 03 ATEX 3565 X

### Introduction

The turbidity sensors METTLER TOLEDO InPro 81XX/82XX \*/\*/\*/ may be used in accordance with EN 60079-14:1997 in the categories 1 to 3 (Zones 0 to 2) and gas groups IIA, IIB und IIC subject to explosion risk through combustible materials in the range of temperature classes T4 to T6.

Radiation flux	≤ 5,19 mW
Pulse energy	≤ 0,519 μJ
Mean radiation flux	≤ 0,13 mW
Mean intensity of irradiation at end of fiber-optic light guide	≤ 0,458 mW/mm2
Irradiation H am at end of fiber-optic light guide	≤ 1,835 μJ/mm2

Turbidity sensor operated in conjunction with transmitter Trb 8300.

### Special conditions

The maximum permissible ambient temperature for each individual temperature class is as shown in the following table:

Temperature class	Max. ambient temperature
T 6	80 °C
T 5	95 °C
T 4	130 °C

- The turbidity sensors METTLER TOLEDO InPro 81XX/82XX \*/\*/\*/ may only be used in the suitable METTLER TOLEDO housings InFit Typ 76\*-\*\*\* or InTrac Typ 7\*\*-\*\*\* in areas subject to explosion risk.
- The measuring system drawings must be followed for installation.
- The metal body of the turbidity sensors as well as the safety weld-in sockets and housings METTLER TOLEDO InFit Typ 76\*-\*\*\* and InTrac Typ 7\*\*-\*\*\* must be electrically connected to the potential equalizing system of the installation.

- The metal body of the turbidity sensors as well as the safety weld-in sockets METTLER TOLEDO housing InFit Typ 76\*\_\*\_\*\_\* and InTrac 7\*\*\_\*\_\*\_\* are, if necessary, to be included into the recurring pressure test of the unit.
- The turbidity sensors may only be operated in conjunction with the transmitter Trb 8300.

### 3. Initial start-up

#### 3.1 Handling

The fiber optic cables have a core of Hard Clad Silica (HCS) fibers. When laying fiber optic cables, the bend radius must always exceed at least 10 cm (4"). Fiberoptic cables may never be shortened independently. Fitting of the SMA plug requires special tooling. The transparent plastic cap on the sensor is to be removed before use. Any formation of deposits of solid particles on the surface is to be avoided. Gas bubbles in the process medium are also able to produce a signal and are therefore to be avoided wherever possible. If this is not viable, please refer to the instruction manual of the METTLER TOLEDO Transmitters Trb 8300 (**chapter: «Calibration»**).

#### 3.2 Installation in vessels/reactors

The turbidity sensors are to be installed at a distance of at least 10 cm (4") from the walls of the vessel, agitator and other internals, in order to prevent disturbances through reflections. When employing special METTLER TOLEDO retractable housings InTrac with modified immersion tubes, the installation guideline valid for each individual housing is to be observed. This step likewise, is to avoid measurement disturbances arising from irrelevant light reflections (**see chapter 7.2: «Accessories/Housings»**).

### 3.3 Installation in pipes

The turbidity sensors are to be installed such that the distance from the pipe walls is at least 10 cm (4"). If this is not feasible in pipes of small diameter, then the sensor can be installed sloped, at an angle of 15...45° or, alternatively, in a knee bend. The pipe must always be completely full of liquid.

### 3.4 Calibration in open laboratory vessels

If the turbidity measuring system is calibrated off-line, it is recommended to use the METTLER TOLEDO calibration tool, «CaliCap» (order no. 52 800 210), for fitting on the sensor shaft. If «CaliCap» is not used, a minimum distance of 10 cm (4") to the vessel walls has to be observed as well. During calibration, the sample should be stirred in order to prevent the sedimentation of solid particles.

### 3.5 Fiber optic cable extensions

Fiber optic cables are available according to the list in **chapter 7.2 «Accessories»**. The duplex cable extensions are supplied complete with four fitted SMA plugs and two couplings. The two SMA plugs of the sensor fixed cable are to be screwed hand-tight onto the couplings. Two SMA plugs of the extension cable can then also be screwed by hand onto the couplings in any order. To protect the cable connection from dust and splashwater, it is recommended to use the METTLER TOLEDO Connection Box IP65 (Nema 4X) (order no. 52 800 241).

If extension cables are installed subsequently at a later date, the measuring system has to be newly calibrated (see instruction manual for the METTLER TOLEDO Transmitters Trb 8300, **chapter: «Calibration»**).

## **4. Maintenance and troubleshooting**

### **4.1 Conditions of warranty**

METTLER TOLEDO guarantees the quality of materials and workmanship within a narrow range of manufacturing tolerances, so that the product purchased is free from any substantial deviations from material and manufacturing quality standards. The warranty is valid for the period of one year from date of delivery ex works. If within this warranty period, any repair or replacement should become necessary, and such cause is not due to misuse or incorrect application, please return the sensor, carriage paid, to your appropriate METTLER TOLEDO agency. Repair work will be carried out free of charge. Final decision on whether the defect is due to a manufacturing error or to incorrect operation of the sensor by the customer is made at the option of the Customer Service department of METTLER TOLEDO. After expiry of the period of warranty, faulty sensors will be repaired or replaced on an exchange basis against payment of the costs involved.

### **4.2 Maintenance**

Dirty or contaminated sensors can deliver incorrect measurement values. If fouling is presumed, the sensor is to be removed and cleaned. Suitable cleaning solutions are mild detergents or strongly diluted acids (< 0.5 % by wt.) such as hydrochloric acid.

When handling acids, precautionary measures are to be taken at all times.

## **5. Removal from operation, storage, disposal**

### **5.1 Removal from operation**

The sensor is only conditionally subject to aging. When employed correctly as intended and appropriately maintained and serviced, the lifetime of the sensor can extend to several years.

Before removing the sensor from the process/process adapter, it must be ensured that the process pressure has been reduced to a safe level and the process temperature lowered to a safe range. Any escape of hot process medium under pressure can cause damage to material/equipment or injury to persons.

After demounting, the sensor should first be flushed with clear water.

If the sensor incurs a defect, it can as a rule not be repaired and must be disposed under observance of prevailing regulations.

### **5.2 Storage**

If the sensor is not needed for use, it can be stored dry.

### **5.3 Disposal**

Disposal is to be carried out by the user in line with valid local regulations. The user must either deliver the sensor to an authorized private or public collection undertaking, or dispose of it himself in accordance with the regulations for such items/materials.

#### **The packaging consists of:**


- cardboard
- foamed plastic

#### **The sensor is comprised of:**


- shaft, O-ring, sealing and window material as per specification
- HCS fiber optic cable, PVC sheathing
- protective tube of DIN 1.4435 (316L) or polyamide

## 6. Dimensional drawings


**InPro® 8200/S(H)**

**InPro® 8100/S**

**InPro® 8050**

## 7. Ordering information

### 7.1 Sensors

Designation	Order no.
InPro8050	52 800 209
InPro8100/S/120	52 800 205
InPro8100/S/205	52 800 206
InPro8100/S/297	52 800 207
InPro8100/S/407	52 800 208
InPro8200/S/Epoxy/120	52 800 216
InPro8200/S/Epoxy/205	52 800 217
InPro8200/S/Epoxy/297	52 800 218
InPro8200/S/Epoxy/407	52 800 219
InPro8200/H/Epoxy/120	52 800 220
InPro8200/H/Epoxy/205	52 800 221
InPro8200/H/Epoxy/297	52 800 222
InPro8200/H/Epoxy/407	52 800 223
InPro8200/S/Kalrez-FDA/120	52 800 224
InPro8200/S/Kalrez-FDA/205	52 800 225
InPro8200/S/Kalrez-FDA/297	52 800 226
InPro8200/S/Kalrez-FDA/407	52 800 227
InPro8200/H/Kalrez-FDA/407*	52 800 215

\* shaft material Hastelloy C22

### 7.2 Accessories

**O-rings, for all sensors:**

Standard O-ring	Additional O-rings	Order no.
Viton®-FDA (10.77 x 2.62 mm)		20 302 1000
	Silicon-FDA (10.77 x 2.62 mm)	20 301 1136
	Kalrez® (10.77 x 2.62 mm)	20 304 1000
	Kalrez®-FDA (10.77 x 2.62 mm)	20 304 1034
	EPDM-FDA (10.77 x 2.62 mm)	20 303 1206

**Fiber optic cable extensions including two couplings, for all sensors:**

<b>Designation</b>	<b>Order no.</b>
Fiber-Kit, 3 m (10 ft)	52 800 228
Fiber-Kit, 5 m (15 ft)	52 800 229
Fiber-Kit, 6 m (18 ft)	52 800 230
Fiber-Kit, 10 m (30 ft)	52 800 231
Fiber-Kit, 15 m (45 ft)	52 800 232
Fiber-Kit, 20 m (60 ft)	52 800 233
Fiber-Kit, 25 m (75 ft)	52 800 234
Fiber-Kit, 30 m (100 ft)	52 800 235
Fiber-Kit, 40 m (130 ft)	52 800 236
Fiber-Kit, 45 m (150 ft)	52 800 172
Fiber-Kit, 50 m (165 ft)	52 800 237
Fiber-Kit, 65 m (200 ft)	52 800 238
Fiber-Kit, 75 m (230 ft)	52 800 177
Fiber-Kit, 90 m (275 ft)	52 800 239
Fiber-Kit, 100 m (330 ft)	52 800 154
Fiber-Kit, 125 m (380 ft)	52 800 158
Fiber-Kit, 170 m (520 ft)	52 800 196

**Additional accessories, for all sensors:**

<b>Designation</b>	<b>Order no.</b>
«CaliCap»	52 800 210
Couplings	52 800 240
Coupling box IP 65	52 800 241
Swagelok adapter NPT 1/2" <sup>1)</sup>	52 800 242

<sup>1)</sup> only to be used with InPro® 8200/S/Epoxy Sensors for pressure range 0...60 bar (0...870 psi)

**Housings for InPro® 8100 and InPro® 8200 (examples):**

<b>Designation</b>	<b>Material</b>	<b>H</b>	<b>a</b>	<b>O-ring</b>	<b>Order no.</b>
InFit 761-25 CP/70/4435/Vi	DIN 1.4435	70	120	Viton®	00 761 3042
InFit 761-25 CP/70/C22/Vi	Hastelloy C22	70	120	Viton®	00 761 3105
InFit 761-25BTB/70	DIN 1.4435	70	120	Silicon-FDA	00 761 3093
InFit 761-25 CIP/70 3.1B	DIN 1.4435	70	120	Silicon-FDA	52 400 491

InFit 761-25 CIP/T 3.1B	DIN 1.4435	–	120	Silicon-FDA	52 400 493
InFit 761-25 CIP/TS 3.1B	DIN 1.4435	–	120	Silicon-FDA	52 400 502
InFit 761-25 CIP/70 TC 1.5"	DIN 1.4435	70	120	Silicon-FDA	52 400 494
InFit 761-25 CIP/70 TC 2"	DIN 1.4435	70	120	Silicon-FDA	52 400 495
InTrac 779-M	DIN 1.4404	70	205	Viton®-FDA	52 403 236
InTrac 799-M	DIN 1.4404	92	297	Viton®-FDA	52 403 237

(H = immersion length)

## Housings for InPro 8050:

Designation	Material	H	a	O-ring	Order no.
InFit 761-25 B PVC	PVC	70	120	Viton®-FDA	52 400 310
InFit 761-25 B PVDF	PVDF	70	120	Viton®-FDA	52 400 311
InFit 761-25 B PP	PP	70	120	Viton®-FDA	52 400 316
InDip 550/1000/PVC/Pg 13.5	PVC	1000	120	Viton®-FDA	52 400 320
InDip 550/1500/PVC/Pg 13.5	PVC	1500	120	Viton®-FDA	52 400 582
InDip 550/2000/PVC/Pg 13.5	PVC	2000	120	Viton®-FDA	52 400 588
InDip 550/2500/PVC/Pg 13.5	PVC	2500	120	Viton®-FDA	52 400 594
InDip 550/3000/PVC/Pg 13.5	PVC	3000	120	Viton®-FDA	52 400 600

(H = immersion length)

## 8. Certificates

### 8.1 Material certificate

Each sensor with a stainless steel shaft made of DIN 1.4435 or 1.4404 (316L) is delivered with a set of 3.1B certificates (complying with EN 10204.3/1.B). All wetted metal parts are identified with an engraved symbol corresponding to the heat number on the paper certificate delivered with the sensor.

Each wetted metal part is polished in order to get a surface roughness lower than 0.4 µm (16 µ in). This represents a roughness grade number of N5 (according to ISO1320:1992).

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[illegible]

[illegible]

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

- BR** **Mettler-Toledo Ind. e Com. Ltda.**, Alameda Araguaia, 451 - Alphaville, BR-06455-000 Barueri/SP, Brazil  
Phone +55 11 4166 74 00, Fax +55 11 4166 74 01
- CH** **Mettler-Toledo (Schweiz) AG**, Im Langacher, CH-8606 Greifensee, Switzerland  
Phone +41 44 944 45 45, Fax +41 44 944 45 10
- D** **Mettler-Toledo GmbH**, Prozeßanalytik, Ockerweg 3, D-35396 Gießen, Germany  
Phone +49 641 507-333, Fax +49 641 507-397
- F** **Mettler-Toledo Analyse Industrielle Sàrl**, 30 Bld. de Douaumont, BP 949, F-75829 Paris Cedex 17, France  
Phone +33 1 47 37 06 00, Fax +33 1 47 37 46 26
- USA** **Mettler-Toledo Ingold Inc.**, 36 Middlesex Turnpike, Bedford MA 01730, USA  
Phone +1 781 301-8800, Fax +1 781 271-0681