Palindrome™

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Sterilized with ethylene oxide. Sterile and non-pyrogenic in unopened and undamaged package.

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

READ ALL INSTRUCTIONS, WARNINGS, AND CAUTIONS CAREFULLY PRIOR TO USE.

DESCRIPTION:

The Palindrome H (heparin coated) 14.5 Fr/Ch cuffed catheter is a radiopaque, urethane catheter with felt cuff, dual extensions and heparin coating.

Each extension has an in-line clamp and luer-lock adapter that is color coded: red for “arterial” outflow of blood, blue for “venous” return. Sealing caps are included with the catheter.

The Palindrome H catheter contains a heparin coating technology applied to the external catheter surfaces (extending from the cuff of the device to the tip) and to the internal catheter surfaces (extending from the luer adapters to the tip).

The coating technology provides a unique biocompatible surface coating involving heparin to reduce platelet adhesion on the catheter.

DIMENSIONS:

CATHETER O. D.	14.5 Fr/Ch (4.8 mm)			
OVERALL LENGTH	36 cm	40 cm	45 cm	50 cm
IMPLANT LENGTH	19 cm	23 cm	28 cm	33 cm

PALINDROME H CATHETER FLOW TABLE WITH SIDE SLOTS

FLOW RATE (ml/min)	Implant Lengths							
	19 cm		23 cm		28 cm		33 cm	
	+VP*	-AP**	+VP*	-AP**	+VP*	-AP**	+VP*	-AP**
200	18	-19	21	-23	23	-24	26	-28
250	24	-28	29	-32	33	-35	37	-40
300	33	-38	39	-44	45	-48	51	-54
350	44	-50	50	-58	59	-64	67	-72
400	56	-64	65	-74	75	-82	83	-93
450	71	-80	79	-92	91	-103	102	-117
500	87	-98	96	-112	108	-126	120	-143

*positive (+VP) Venous Pressure (mm Hg)
**negative (-AP) Arterial Pressure (mm Hg)

INDICATIONS:

The Palindrome H 14.5 Fr/Ch Chronic Hemodialysis Catheter with heparin coating is intended for hemodialysis, apheresis, and infusion.

The advantage of the heparin coating on this catheter in reducing platelet adhesion to the catheter surface was supported by bench and animal testing. It may be inserted either percutaneously or by cutdown.

TEST SUMMARY:


The advantage of the heparin coating on this catheter in reducing platelet adhesion to the catheter surface is supported by:

- A two hour circulating blood loop test demonstrating a 60% reduction in platelet adhesion on the catheter surface at p<0.05.
- A coating durability test, where the catheter was subjected to 720 hours of simulated dialysis conditions, and maintained heparin activity levels at twice the minimum activity level required to achieve a 60% reduction in platelet adhesion.
- An in-vivo ovine model using six sheep (periodically perfused to simulate dialysis for 24 days) where the reduction in thrombus formation was 82% at p<0.05. Disseminated fibrin sheath and clotting was more commonly observed along the shaft of the uncoated catheter than the coated catheter.

TEST METHOD DETAILS:


In-vitro evaluations of the coated catheters were performed using a test model which incorporates fresh heparinized bovine blood to assess the relative thromboresistance of the coated catheter as compared to a non-coated catheter. The blood, with radiolabeled autologous platelets, was circulated for 2 hours.

Retrieved catheters were visually inspected and then placed in a gamma counter for quantification of platelet adhesion on the catheter surface. The radioactivity data demonstrates that the coated catheter had 60% less platelets adhered to the surface compared with the uncoated catheter.


Total end-point platelet accumulation normalized to the uncoated control.

In-vivo evaluations of the coated catheters were performed using an ovine model. The testing was conducted on 6 anti-coagulated sheep with a coated and non-coated catheter implanted into the right and left jugular veins of the same sheep. Routine weekly blood perfusion sessions were performed on both catheters to simulate dialysis. Gravimetric analysis performed on the thrombus extracted from the external surfaces of both the coated and non-coated catheters demonstrated an 82% reduction in total thrombus formation on coated catheters after an average of 24 days of implantation as compared to non-coated catheters. Disseminated fibrin sheath and clotting was more commonly observed along the shaft of the uncoated catheter than the coated catheter.

The durability of the coating was assessed in an in-vitro test model that simulates the dynamic flow environment of a dialysis session. The model involves 37°C Saline flowing through the internal surfaces and around the external surfaces of the catheter for a time period that simulates over 12 months of dialysis sessions on the ID of the catheter and over 30 days on the OD of the catheter. The chart below shows that between 60% and 70% of the Palindrome H catheter heparin activity remains after 720 hours of continuous flow. This heparin activity is significantly above the minimum heparin activity established during in-vitro blood flow evaluations to achieve a 60% reduction in platelet adhesion.

* The established efficacy threshold was determined in an in-vitro circulating bovine blood model using coated catheters with varying levels of heparin activity. The blood, with radiolabeled autologous platelets, was circulated for 2 hours. Platelet counts were quantified for each of the coated catheters with varying heparin activity levels and compared to the uncoated catheter. The results demonstrated that a catheter with 43% of the Palindrome H catheter heparin activity still provides a 60% reduction in platelet adhesion on the catheter surface.

CONTRAINDICATIONS:

Do not use this catheter in thrombosed vessels or for subclavian puncture when ventilator is in use.

Heparin coated catheters should not be used in individuals with documented hypersensitivity to heparin or porcine based products. Heparin coated catheters should not be used in patients with severe thrombocytopenia, uncontrollable active bleeding disorders, or with skin necrosis from previous heparin use.

POTENTIAL COMPLICATIONS:

Potential complications include:

sepsis; thrombosis/stenosis of vein; exit site infection; cardiac arrhythmia; air embolism; subcutaneous tunnel infection; hemorrhage; hemothorax; pneumothorax; hematoma; cardiac tamponade; trauma to major vessel or right atrium; brachial plexus injury; catheter thrombosis; retroperitoneal bleed; femoral nerve damage; femoral artery damage; femoral artery dissection; femoral vein occlusion; lower extremity ischemia; pulmonary embolism; pleural injury; arterial puncture; mediastinal widening; deep vein thrombosis of the lower extremity; hemomediastinum; recurrent laryngeal nerve palsy; dissection or occlusion of the carotid artery; thoracic duct laceration; hypersensitivity reactions; local irritation; skin necrosis; endocarditis.

WARNINGS AND PRECAUTIONS:

- It is important to consult the Summary of Product Characteristics for heparin sodium before using this product.
- Although the Palindrome H catheter incorporates a firmly bonded coating, care should be exercised to avoid excessive abrasion of the surface.
- The heparin coating used on the Palindrome H catheter is not intended as a replacement for priming the catheter lumens with heparinized saline. Standard priming procedures still apply.
- As with any heparin-based product, the following conditions should be considered:
 - Hypersensitivity: Patients with documented hypersensitivity to heparin or porcine based products should not receive heparin coated catheters.
 - Thrombocytopenia: Heparin-induced thrombocytopenia has been reported with the use of heparin-coated catheters.
- The use of 46.7% citrate as a locking solution is not recommended. If used as a locking solution it may compromise the efficacy of the internal lumen coating activity.
- The catheter should be inserted and removed only by a qualified, licensed physician or other healthcare practitioner authorized by and under the direction of such physician.
- The medical techniques and procedures described in these instructions do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient.
- Observe sterile technique at all times when handling or using the catheter.
- Do not use the catheter if package has been previously opened or is damaged.
- Do not use the catheter if it appears damaged or defective.
- Avoid air embolism by keeping catheter extension tubing clamped at all times when not in use and by filling the catheter with sterile saline prior to implantation. With each tubing change, purge air from the tubing and aspirate any air in catheter.
- Do not clamp extension tubing with insertion stylets loaded as it may result in stylet damage.
- Percutaneous insertion of the catheter in the subclavian vein can be technically difficult. The right internal jugular is preferable.
- To avoid vessel perforation and damage, do not insert the guidewire, dilators, or valved pull-apart sheath/introducer forcibly.
- Do not insert the valved pull-apart sheath/introducer further than necessary: depending upon patient size and access site, it may not be necessary to insert the entire length of the introducer into the vessel.
- The valved pull-apart sheath/introducer is designed to reduce blood loss and the risk of air intake but it is not a hemostasis valve.

- The valved pull-apart sheath/introducer is not intended to create a complete two-way seal nor is it intended for arterial use.
 - The valve will substantially reduce air intake. At -12 mm Hg vacuum pressure the valved pull-apart sheath/introducer may allow up to 4 cc/sec of air to pass through the valve.
 - The valve will substantially reduce the rate of blood flow but some blood loss through the valve may occur.
- Use the guidewire straightener to insert the “J” end of the guidewire into the introducer needle. Do not insert or withdraw the guidewire forcibly from any component; the wire could break or unravel.
- Do not nick the catheter when suturing.
- Do not tie the suture too tightly at the venotomy site.
- Prolonged exposure to ultraviolet light can damage the catheter.
- Do not use acetone on any part of the catheter. Aqueous-based povidone iodine, Exsept™, Hibiclens™ (Chlorohexidine), amikun 50%, hydrogen peroxide, Neosporin™ antibiotic ointment, bacitracin ointment, Bactroban™ cream, isopropyl alcohol 70%, Chloraprep™ can be used. Inter-mixing of these solutions has not been tested and is not recommended.
- Overtightening catheter connections can crack some adapters.
- Do not clamp the dual lumen portion of the catheter; clamp only the extensions. Use only smooth-jawed forceps for clamping when not using the clamp supplied with the catheter.
- Clamping the catheter repeatedly in the same spot could weaken the tubing: change the position of the clamp regularly to prolong the life of the tubing. Avoid clamping near the adapter and hub.
- Exercise caution when using sharp instruments near the catheter. Catheter tubing can tear when subjected to nicks, excessive force, or rough edges.
- Inspect the catheter frequently for nicks, scrapes, cuts, etc. which could impair its performance.
- When infusing heparin, flush quickly and clamp immediately to ensure heparin reaches the distal end of the lumen. Do not infuse against a closed clamp or forcibly infuse a blocked catheter: back pressure could force the adapter out of the tubing.
- Remove the catheter as soon as it is no longer necessary.
- When removing the catheter, DO NOT use a sharp, jerking motion or undue force; this may tear the catheter. Free the cuff and surfaces from the tissue prior to removal.
- Discard the catheter after single use. Do not resterilize.
- This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

PRECAUTIONS:

- Consult pharmacy or appropriate literature for compatibility data before using heparin coated catheters with any drug product. Drugs that bind or interact with heparin should not be infused through the lumens of the heparin coated catheter.
- There is increased risk of bleeding with anticoagulants, platelet inhibitors or any other drug which may interfere with coagulation, including some cephalosporins. Drugs that interfere with platelet aggregation should be used with care. There is an increased risk of haemorrhage with ketorolac. Reduced activity of heparin has been reported with simultaneous intravenous glycerol trinitrate infusion.
- Heparin may also interfere with diagnostic tests. Pseudo-hypocalcaemia (in haemodialysis patients), artifactual increases in total thyroxine and triiodothyronine, simulated metabolic acidosis and inhibition of the chromogenic lysate assay for endotoxin may all result from interaction with heparin. Immunoassay determination of aminoglycosides may also be affected.

RECOMMENDED

- Use only luer-lock (threaded) connections (including syringes, bloodlines, and sealing caps) with the catheter's adapters.
- Use a straight tip guidewire when removing a Palindrome H catheter with side slots for an over-the-guidewire catheter exchange.

INSERTING THE PALINDROME H 14.5 Fr/Ch CUFFED CATHETER

Sterile Components Supplied in Kit

Catheter, Tal VenaTrac Insertion Stylets (optional), 18 gauge (1.2 mm) J Introducer Needle, Tunneler, 12 mL Syringe, #11 Scalpel, 0.038 in (0.965 mm) J/Straight Guidewire, 4" x 4" Gauze Sponges, 14 Fr (4.7 mm) Dilator, Wound Dressing, 12 Fr (4.0 mm) Dilator, Sealing Caps, 16 Fr (5.3 mm) Valved pull-apart Sheath/Introducer.

Sterile Components Supplied in Sport Pack

Catheter, Insertion Stylets (optional), 14 Fr (4.7 mm) Dilator, Tunneler, 12 Fr (4.0 mm) Dilator, Sealing Caps, 16 Fr (5.3 mm) Valved pull-apart Sheath/Introducer.

Please reference product label for exact components supplied.

Additional Sterile Components that may be Required

Prepping Agents, Normal Saline, Drapes, Needle Holder, Mask, Gloves, Gown, Suture with Curved Needle, Syringes/Needles, Razor, Local Anesthetic, Heparin (in concentrations approved by your institution).

NOTE: Components needed may vary depending on insertion procedure used. The Palindrome H 14.5 Fr/Ch cuffed catheter may be inserted either percutaneously or by cutdown.

INSERTION SITE

The Palindrome H (heparin coated) 14.5 Fr/Ch cuffed catheter is ideally placed in the right atrium via the right internal jugular. While the catheter also can be placed in the external jugular, subclavian, or saphenous vein, the right internal jugular is strongly recommended for the following reasons:

- The internal jugular permits easier positioning of the catheter tip in the right atrium.
- The size and location of the external jugular makes insertion difficult.
- Use of the subclavian vein for catheter placement may result in subclavian vein stenosis. Subclavian vein stenosis may prevent the future use of the ipsilateral extremity for permanent access.
- Subclavian vein placements are at risk for higher insertion complications
- The saphenous vein should be used only as a last alternative due to the possibility of insertion complications.

For optimal catheter function: Per Kidney Dialysis Outcome Quality Initiatives (K/DOQI), the catheter tip must be adjusted to the level of the caval atrial junction or beyond to ensure optimal blood flow. For this reason, it is preferable to insert the catheter on the patient's right side. For large patients, and for patients whose right-side veins are unusable, 40 cm, 45 cm, 50 cm catheters are available for placement via the left jugular or subclavian vein.

PREPARATION

The operating room or IR suite are the preferred locations for catheter placement. **Both cutdown and percutaneous procedures require confirmation of correct placement by fluoroscopy or chest x-ray.**

1. Provide a sterile operative field: use sterile drapes, instruments, and accessories. Perform surgical scrub. Wear a gown, cap, gloves, and mask. Have the patient wear a mask.
2. Place the patient in a supine position and expose the upper chest or the groin side to be accessed.

For subclavian and jugular insertion: Turn the patient's head slightly to the side to expose the insertion site. The Trendelenberg position may facilitate insertion.

3. Shave the access site (optional) and prep the area in the established manner. Isolate the access site with sterile drapes.

NOTE: When performing a subclavian insertion on a patient with large breasts, it is best to draw landmarks while the patient is sitting up to prevent catheter tip migration.

4. Fill the catheter with sterile heparinized saline and clamp the extensions immediately.
5. Administer local anesthetic to the skin and underlying tissue at the insertion site.

WARNING: To prevent air embolism, keep the catheter clamped at all times when not attached to a syringe, IV tubing, or bloodlines.

NOTE: If using insertion stylets, catheter will not be filled until after placement. Load stylets into catheter lumens such that a minimum length of 8 cm extends distal from the catheter adapters, in order to allow for proper catheter tunneling.

PERCUTANEOUS PROCEDURE MODIFIED SELDINGER (LITTLEFORD-SPECTOR) TECHNIQUE

Cannulating the Vessel


1. Flush an 18 gauge introducer needle with heparinized normal saline. Insert the needle through the primary incision and advance it into the vein, in the direction of blood flow. Aspirate a small amount of blood to ensure the needle is correctly positioned in the vein.

CAUTION: If arterial blood is aspirated, remove the needle and apply immediate pressure to the site for at least 15 minutes. Ensure that the bleeding has stopped and that no hematomas have developed before attempting to cannulate the vein again.

- 2a. Disconnect the syringe from the needle and promptly insert the flexible “J” end of the guidewire through the introducer needle. Failure to insert the wire promptly may lead to blood loss through the needle. Advance the wire into the vein.

CAUTION FOR JUGULAR AND SUBCLAVIAN INSERTION: The length of wire inserted is determined by the size of the patient. Cardiac arrhythmia may result if the guidewire passes into the right atrium. If symptoms occur, pull back the guidewire until they disappear. If the guidewire meets resistance, do not pull it back through the needle. Remove the wire and the needle together as a unit and begin again with new needle and guidewire (Figure 1).

NOTE: If utilizing the stylets for insertion, the guidewire provided is recommended. Otherwise, greater than a 0.035 inch hydrophilic or greater than a 0.038 inch stainless steel wire is contraindicated.


- 2b. Withdraw the introducer needle, leaving the guidewire in the vein.

CREATING THE SUBCUTANEOUS TUNNEL

NOTE: For ease in dressing the exit site and for patient comfort, locate the subcutaneous tunnel below the entry site. A tunnel with a wide, gentle arc lessens the risk of kinking at the cuff. The tunnel should be short enough to keep the Y-hub of the catheter from entering the exit site, yet long enough to keep the cuff a minimum of 2 cm from the exit site.

1. Make a small incision at the insertion site. (The right internal jugular is the preferable site for percutaneous placement as subclavian placement is technically difficult.) Make a second incision parallel to the first at the exit site. Make the exit incision just long enough to accommodate the cuff, approximately 1 cm.
2. Use blunt dissection to create the subcutaneous tunnel.
 - a. Attach the catheter to the tunneler by sliding the catheter tip onto the bifurcated tunneler tines until the catheter tip meets the base of the tines. Slide the sheath completely over the connection until it stops, being careful that the sheath smoothly transitions over the catheter tip (see Figure 2). If desired, bend the tunneler into a wide arc to make a curved tunnel.

- b. Insert the tunneler into the exit site and create a short subcutaneous tunnel, emerging at the insertion site. The catheter will thread through the tissue as the tunnel is created. Ensure the catheter passes through the tunnel to the primary insertion site.
- c. Remove the catheter carefully from the tines. To do so, slide the sheath back, grasp the tip of the catheter, and gently pull the catheter from the tunneler tines. **Discard the tunneler.**

INSERTING THE CATHETER USING THE VALVED PULL-APART SHEATH/INTRODUCER


OPTIONAL: To ease insertion of the valved pull-apart sheath/introducer, the vein can be pre-dilated with the dilator(s) provided.

1. Thread the dilator(s) over the end of the guidewire and advance it into the vein using a rotating motion to assist passage through the tissue.



CAUTION: Do not force the dilator(s). Ensure that the guidewire does not advance further into the vein.

Remove and discard the dilator(s).

2. Remove the valved pull-apart sheath/introducer assembly from the package; remove dilator from the sheath (Figure 3).

3. Slide the valve over the sheath opening, and insert the dilator through the valve and lock in place using the rotating collar. (Refer to Figures 4 and 5.)

4. a. Thread the locked valved pull-apart sheath/introducer assembly over the end of the guidewire.

CAUTION: To avoid damaging the tissue and the sheath tip, do not let the sheath advance over the dilator. The two must be grasped as one unit.

- b. With a rotational motion, advance the assembly into the vein only as far as necessary.

Do not force the introducer into the vessel. Do not insert it further than necessary for the patient's size and access site. Ensure that the guidewire does not move further into the vein.

- c. Hold the sheath in place, unlock the rotating collar, and gently remove the dilator and guidewire together (discard the dilator and the guidewire).
5. Advance catheter through the valve. To prevent kinking the catheter, it may be necessary to advance in small steps by grasping the catheter close to the sheath (Figure 6).

WARNING: Ensure the catheter is filled with heparinized saline and is free of air bubbles before inserting it into the vein.


Figure 6

6. Continue inserting the catheter through the valved sheath and into the vein. For optimal catheter function: Per Kidney Dialysis Outcome Quality Initiatives (K/DOQI), the catheter tip must be adjusted to the level of the caval atrial junction or beyond to ensure optimal blood flow (Figure 8).
7. Aspirate to verify patency and clamp the extension.
8. After the catheter has been positioned, grasp both tabs firmly and crack the sheath handle in half.

NOTE: Steps 9 and 10 refer to Figure 7.

9. Peel the non-valved side of the handle partially away from the catheter.
10. Near the valve, hold the catheter firmly in position and pull the valve off of the catheter.

NOTE: Some resistance will be experienced while pulling the valve off of the catheter.


Figure 7

11. Remove the sheath from the patient by holding the catheter in place and pulling the separated tabs away from the entry site simultaneously at a 180° angle. The sheath will separate from the catheter, leaving the catheter in the vein.

CAUTION: Do not allow the catheter to move out of the vein with the sheath. Ensure that the vein is not bleeding around the catheter.

12. Use fluoroscopy or portable x-ray to view the catheter. Position the catheter correctly in the right atrium.


Figure 8: Implanted Catheter

13. Confirm correct placement and catheter function by aspirating venous blood from both the arterial and venous lumina, then flush 5 mL of sterile normal saline into each lumen. Follow with heparinized saline. Clamp the extensions immediately utilizing the positive pressure technique. (See priming volumes under HEPARINIZATION.) Attach a sterile sealing cap to each adapter.

WARNING: To prevent air embolism, expel all air from syringes before injecting solutions.

14. Suture the entry site. The exit site should not require suture.
15. Suture the catheter hub to the skin utilizing a prolene or nylon 3-0 or 4-0 suture on the extension wings.

CAUTION: Do not suture through any part of the catheter. Remove skin suture by the 4th or 5th day to diminish erythema in the area. The catheter can also be immobilized with either gauze or transparent dressings..

16. Apply dressing to catheter exit site and insertion site.

SHEATHLESS INSERTION OR CATHETER EXCHANGE USING TAL VENATrac™ OVER-THE-WIRE INSERTION STYLETS

(Neither the Palindrome H Catheter nor the Tal VenaTrac insertion stylet are intended for femoral placement)

1. Load stylets fully into catheter and lock stylet hubs to lumen adapters. Ensure that the white flat side of the longer stylet containing the stylet sidehole is oriented against the septum of the catheter lumen (Figure 9).

NOTE: Flush the Tal VenaTrac insertion stylets with sterile saline prior to use. Wetting the outer surfaces of the stylets as well may facilitate advancement into the catheter.


Figure 9

OPTIONAL: To ease insertion of the catheter, the vein can be pre-dilated with the dilator(s) provided.

2. Thread the dilator(s) over the end of the guidewire and advance it into the vein using a rotating motion to assist passage through the tissue.

CAUTION: Do not force the dilator(s). Ensure that the guidewire moves freely within the dilator(s) and does not advance further into the vein.


Remove and discard the dilator(s).

3. Thread the guidewire 1 cm into the distal tip of the longer stylet such that it exits the stylet sidehole (Figure 9). To assist the guidewire exit, the longer stylet tip should be flexed 45°.
4. Continue threading the wire into the distal tip of the adjacent shorter stylet (Figure 9) and advance until the wire passes out the proximal end of the catheter-stylet assembly.

- While maintaining control of the guidewire, pinch the catheter-stylet slightly, approximately 2 cm from the distal tip, and advance into the vein (Figure 10).

NOTE: For difficult placements, if so desired, the catheter may be rotated in a clockwise direction to ease advancement. If excessive resistance is felt, stop advancement and retract the entire assembly. Verify that the guidewire has not kinked and re-attempt catheter advancement.

CAUTION: Do not force the assembly. Ensure that the guidewire does not advance further into the vein.

- Continue to pinch slightly and advance the catheter-stylet assembly until catheter tip is positioned as desired (Figure 10). Per Kidney Dialysis Outcome Quality Initiatives (K/DOQI), the catheter tip must be adjusted to the level of the caval atrial junction or beyond to ensure optimal blood flow (Figure 8).
- Holding the catheter in place, gently remove the guidewire and discard. Next, immediately unlock, remove and discard each of the two Tal VenaTrac insertion stylets.

CAUTION: Do not remove insertion stylets while guidewire is in place as it may result in catheter damage. Ensure that the catheter does not move out of the vein while removing insertion stylets.

CAUTION: If guidewire meets resistance, do not pull it back through the catheter-stylet assembly. Remove the catheter, insertion stylets, and wire together as a unit and begin again with a new catheter.

- Aspirate to verify patency and clamp the extension.
- Use fluoroscopy or portable x-ray to view the catheter. Position the catheter correctly in the right atrium (Figure 8).
- Confirm correct placement and catheter function by aspirating venous blood from both the arterial and venous lumina, then flush 5 mL of sterile normal saline into each lumen. Follow with heparinized saline. Clamp the extensions immediately utilizing the positive pressure technique. (See priming volumes under HEPARINIZATION.) Attach a sterile sealing cap to each adapter.

WARNING: To prevent air embolism, expel all air from syringes before injecting solutions.

- Suture the entry site. The exit site should not require suture.
- Suture the catheter hub to the skin utilizing a prolene or nylon 3-0 or 4-0 suture on the extension wings.

CAUTION: Do not suture through any part of the catheter. Remove skin suture by the 4th or 5th day to diminish erythema in the area. The catheter can also be immobilized with either gauze or transparent dressings.

- Apply dressing to catheter exit site and insertion site.

CUTDOWN PROCEDURE

- Make an adequate incision following the skin lines over the desired vessel.
- Ensure the chosen vein is large enough to accept the catheter.
- If using a large vein such as the jugular, immobilize the vein and place a purse string suture around the venotomy site. If using a smaller vein such as the saphenous, tie down the distal ligature, but leave the proximal one untied until the catheter is in position.

- Make a counter incision for the exit site of the subcutaneous tunnel. Make the incision just long enough to accommodate the cuff, approximately 1.0 cm.

NOTE: For ease in dressing the exit site and for patient comfort, locate the subcutaneous tunnel below the entry site. A tunnel with a wide, gentle arc lessens the risk of kinking the catheter at the cuff. The tunnel should be short enough to keep the Y-hub of the catheter from entering the exit site, yet long enough to keep the cuff a minimum of 2cm from the exit site.

- Use blunt dissection to create a short subcutaneous tunnel, emerging at the cutdown incision.

- Attach the catheter to the tunneler by sliding the catheter tip into the bifurcated tunneler tines until the catheter tip meets the base of the tines. Slide the sheath completely over the connection until it stops, being careful that the sheath smoothly transitions over the catheter tip (see Figure 2 in Percutaneous Procedure section). If desired, bend the tunneler into a wide arc to make a curved tunnel.
 - Insert the tunneler into the exit site and create a short subcutaneous tunnel, emerging at the cutdown incision. The catheter will thread through the tissue as the tunnel is created. Ensure the catheter passes through the tunnel to the cutdown incision.
 - Remove the catheter carefully from the tines. To do so, slide the sheath back, grasp the tip of the catheter, and gently pull the catheter from the tunneler tines. **Do not tug the catheter tubing.** Discard the tunneler.
- Use forceps to grasp the designated venotomy vein transversely across its entire width.
 - Insert a #11 blade at the midpoint of the vein diameter; cut anteriorly for a controlled venotomy. When the vein is released from the forceps, it will display an oval venotomy. This may be dilated further using a mosquito hemostat.

WARNING: Ensure the catheter is filled with heparinized saline and is free of air bubbles before inserting it into the vein.

- Grasp the end of the catheter with forceps and insert it into the vein. Slacken the proximal ligature to allow the catheter to pass with minimal back bleeding. Advance the catheter into the vein: for subclavian and jugular insertion, advance the catheter tip into the right atrium. For optimal catheter function: Per Kidney Dialysis Outcome Quality Initiatives (K/DOQI), the catheter tip must be adjusted to the level of the caval atrial junction or beyond to ensure optimal blood flow. (See Figure 8 in Percutaneous Procedure section.)
- Pull the purse string suture (or proximal ligature) closed, but do not tie it before determining the exact position of the catheter.
- Use fluoroscopy or portable x-ray to view the catheter. Move or reposition the catheter until it is positioned correctly in the right atrium as indicated in Figure 8 in Percutaneous Procedure section.
- Tie the purse string suture (or proximal ligature) snugly around the catheter. Tie it just tightly enough to control bleeding at the venotomy; do not occlude the catheter.

CAUTION: Do not allow the cuff to enter the vein or the venotomy. Ideally, it should not touch the vein.

- Confirm correct placement and catheter function by aspirating venous blood from both the arterial and the venous lumina, then flush 5 mL sterile normal saline into each lumen. Follow with heparinized saline. Clamp the extensions immediately utilizing the positive pressure technique. (See priming volumes under HEPARINIZATION.) Attach a sterile sealing cap to each adapter.

WARNING: To prevent air embolism, expel all air from syringes before injecting solutions.

- Suture the entry site. The exit site should not require suture.
- If a skin suture is desired, suture the catheter hub to the skin utilizing a prolene or nylon 3-0 or 4-0 suture on the extension wings

CAUTION: Do not suture through any part of the catheter. Remove skin suture by the 4th or 5th day to diminish erythema in the area. The catheter can also be immobilized with either gauze or transparent dressings.

- Apply dressing to the catheter exit site and cutdown incision.

REMOVING THE PALINDROME H 14.5 Fr/Ch CUFFED CATHETER

To remove the catheter, free the cuff from the tissue and pull the catheter gently and smoothly. Do not use sharp, jerking motions or undue force. This could tear the catheter.

WARNING: If the catheter offers resistance, do not pull further. Perform a cutdown and remove all sutures at the venotomy site.

DISPOSAL

After use, the catheter and accessories are considered a biohazard. Handle and dispose of these in accordance with accepted medical practice and all applicable laws and regulations.

HEPARINIZATION

To maintain patency between dialysis or apheresis treatments, keep the lumina of the catheter filled with the appropriate concentration and volume of heparin. In most cases, 5,000 units/mL is most successful (refer to the priming volume of the catheter). Approved heparin concentrations vary with each institution. Be sure to use those concentrations approved by your facility.

PRIMING VOLUMES

Catheter	Overall Length	Arterial	Venous
ADULT	36 cm	1.6 mL	1.6 mL
	40 cm	1.9 mL	1.9 mL
	45 cm	2.1 mL	2.1 mL
	50 cm	2.3 mL	2.3 mL

Heparinize only after use. Before initiating treatment, aspirate indwelling heparin and discard. After treatment, flush well and instill fresh heparin. If the interdialytic period is less than two days, or if apheresis is performed daily, a lower concentration of heparin may be desirable.

In all cases, the patient’s condition must be considered when choosing a heparin regime. Use less heparin in children and in adults with bleeding disorders.

SUPPLIES

- 10-20 mL syringes
- 3 mL syringe
- 20 gauge (0.91 mm) 2.5 cm needles
- vial heparin (in concentration approved by your institution)
- vial sterile normal saline
- povidone-iodine swabs

PREPARATION

1. Prepare supplies on a clean surface.
2. Wash hands thoroughly with soap and water.
3. Scrub the area surrounding the cap and catheter for 5 minutes with a povidone-iodine swab. **Allow to air dry.**
4. Open syringe and needle packages. Place the needle on the sterile syringe, using aseptic technique.
5. Remove the tops of the saline and heparin vials and swab the injection area with povidone-iodine. **Allow to air dry.**
6. Prepare the appropriate dilute heparin solution.

PROCEDURE

1. Remove sealing cap and aspirate indwelling heparin from the catheter before infusing fresh heparin or initiating treatment.
2. Flush the lumen with 10 to 20 mL sterile normal saline.
CAUTION: Before flushing, pull the plunger back to verify blood flow and to ensure there are no blood clots. Do not flush clots through the catheter (see “Thrombus Formation”).
3. Infuse fresh heparin, flushing quickly to ensure that heparin reaches the distal end of the lumen, and clamp immediately. Infusing or clamping too slowly may cause heparin to exit the catheter from the catheter slot, leaving the distal catheter tip unprotected from thrombus formation. Do not infuse against a closed clamp or forcibly infuse a blocked catheter: back pressure could force the adapter to loosen and potentially come out of the tubing. Perform procedure for both lumina.
Once the lumen has been primed, keep the extension clamped when not attached to a bloodline or syringe. If the extension is unclamped, the priming volume will increase slightly as a result of the tube returning to its “normal” unclamped state. This creates a vacuum at the tip, causing blood to be drawn into the distal portion of the catheter, ultimately resulting in a thrombus.

MANAGEMENT OF ONE-WAY OBSTRUCTION

One-way obstruction, which exists when a lumen can be flushed easily, but blood cannot be aspirated, usually is caused by tip malposition. One of the following adjustments may resolve the obstruction:

- Reposition the patient.
- Have the patient cough.
- Provided there is no resistance, flush the catheter vigorously with sterile normal saline to try to move the tip away from the vessel wall.

The following procedures may require a physician’s order:

- Consider using a thrombolytic agent.
- If the one-way obstruction exists in the arterial lumen, consider reversing the bloodlines. The patient may be dialyzed by connecting the arterial bloodline to the venous adapter and the venous bloodline to the arterial adapter.

THROMBUS FORMATION

NEVER FORCIBLY FLUSH AN OBSTRUCTED LUMEN.
If either lumen develops a thrombus, first attempt to aspirate the clot with a syringe. If aspiration fails, the physician may attempt to lyse the clot with a thrombolytic agent.

WARNING: Thrombolytic agents may cause systemic fibrinolysis if infused into the circulation. Refer to manufacturer’s instructions, indications for use, and contraindications before using thrombolytic agents. Streptokinase is not recommended; it is reported to be anaphylactogenic.

DIALYSIS

For priming volumes, see Heparinization section of these instructions.

1. Prepare dialysis machine and prime bloodlines in the usual manner. Maintain sterile technique when handling or using the catheter. Scrub the adapters, sealing caps, clamps, extension tubes, and Y-connector of the catheter with an aqueous-based povidone-iodine solution.
2. a. Verify that the arterial extension is clamped before removing the sealing cap.
b. Remove the sealing cap from the adapter and attach a luer-lock syringe.
c. Confirm patency of the lumen by aspirating indwelling heparin until venous blood appears. Leave the syringe in place.
d. Close the clamp on the extension.
3. Repeat Steps 2a. through 2d. for the venous lumen.

WARNING: When connecting bloodlines to the catheter, do not permit air to enter the blood path.

4. Remove the syringe and connect the arterial bloodline to the arterial (red) adapter. Open the clamps on the arterial extension and arterial and venous bloodlines and turn on the blood pump.
5. Prime the extracorporeal circuit with the patient’s blood and turn off the blood pump. Ensure the venous extension is clamped, then remove the syringe and connect the venous bloodline to the catheter’s venous (blue) adapter. Open the clamps on the venous extension and bloodline and turn on the pump.
6. Begin treatment.

POST DIALYSIS

Prepare syringes with sterile normal saline and heparin.

1. Stop the blood pump. Close the clamp on the arterial extension and clamp the arterial bloodline at the connection site. Disconnect the arterial bloodline from the adapter of the catheter.
2. Connect a 10-20 mL syringe filled with sterile normal saline to the arterial adapter; open the clamp on the arterial extension and flush the blood from the arterial lumen of the catheter. Reclamp the extension, Heparinize the lumen with the appropriate volume/concentration of heparin.
3. Rinse back the blood in the extracorporeal circuit via the catheter’s venous lumen.
4. After the patient’s blood has been rinsed back, turn off the blood pump. Clamp the venous extension and disconnect the venous bloodline from the venous adapter of the catheter.

5. Connect a 10-20 mL syringe filled with sterile normal saline to the venous adapter. Open the clamp on the venous extension and flush all remaining blood from the venous lumen of the catheter. Reclamp. Heparinize the lumen with the appropriate volume/concentration of heparin.
6. Ensure the clamps are closed on both extensions. Remove syringes and attach a sealing cap to each adapter.

WARNING: Keep the catheter clamped at all times except when connected to bloodlines or syringe during treatment.

CATHETER CARE GUIDES

For further information and a copy of Catheter Care guidelines for the clinician, contact your Covidien representative. In the United States call 1-800-962-9888 and for international inquiries, call 508-261-8000.

ADDITIONAL READING

Moss, Alvin H. et al. "Use of a Silicone Catheter with a Dacron Cuff for Dialysis Short-Term Vascular Access." American Journal of Kidney Diseases, Vol. XII, No. 6 (December), 1988: pp. 492-498.

Schwab, Steve J. et al. "Prospective Evaluation of a Dacron Cuffed Hemodialysis Catheter for Prolonged Use." American Journal of Kidney Diseases, Vol. XI, No. 2 (February), 1988: pp. 166-169

Kirkpatrick, W. G., Culpepper, R. M. & Sirmon, M.D., Frequency of Complications with Prolonged Femoral Vein Catheterization for Hemodialysis Access. Nephron: 1996; 73: pp. 58-62.

Zaleski, G. X., Lorenz, J. M., Garofalo, R. S., Moscate, M. A., Rosenblum, J. D. and Leef, J. A. (1998). Experience with Tunneled Femoral Hemodialysis Catheters. American Journal of Radiology, 172: pp. 493-496.

Heparin Sodium Summary of Product Characteristics
K/DOQI Guidelines, 2007.

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LECTURAS ADICIONALES

Moss, Alvin H. et al. "Use of a Silicone Catheter with a Dacron Cuff for Dialysis Short-Term Vascular Access." *Ó American Journal of Kidney Diseases*, Vol. XII, No. 6 (diciembre), 1988: pp. 492-498.

Schwab, Steve J. et al. "Prospective Evaluation of a Dacron Cuffed Hemodialysis Catheter for Prolonged Use." *Ó American Journal of Kidney Diseases*, Vol. XI, No. 2 (febrero), 1988: pp. 166-169

Kirkpatrick, W. G., Culpepper, R. M. & Sirmon, M.D., Frequency of Complications with Prolonged Femoral Vein Catheterization for Hemodialysis Access. *Nephron*: 1996; 73: pp. 58-62.

Zaleski, G. X., Lorenz, J. M., Garofalo, R. S., Moscate, M. A., Rosenblum, J. D. and Leef, J. A. (1998). Experience with Tunneled Femoral Hemodialysis Catheters. *American Journal of Radiology*, 172: pp. 493-496.

Heparin Sodium Summary of Product Characteristics K/DOQI Guidelines, 2007.

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Moss, Alvin H. et al. "Use of a Silicone Catheter with a Dacron Cuff for Dialysis Short-Term Vascular Access." American Journal of Kidney Diseases, Vol. XII, No. 6 (December), 1988: sid. 492-498.

Schwab, Steve J. et al. "Prospective Evaluation of a Dacron Cuffed Hemodialysis Catheter for Prolonged Use." American Journal of Kidney Diseases, Vol. XII, No. 2 (December), 1988: sid. 166-169

Kirkpatrick, W. G., Culpepper, R. M. & Sirmon, M.D., Frequency of Complications with Prolonged Femoral Vein Catheterization for Hemodialysis Access. Nephron: 1996; 73: sid. 58-62.

Zaleski, G. X., Lorenz, J. M., Garofalo, R. S., Moscate, M. A., Rosenblum, J. D. and Leef, J. A. (1998). Experience with Tunneled Femoral Hemodialysis Catheters. American Journal of Radiology, 172: sid. 493-496.

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Moss, Alvin H. et al. "Use of a Silicone Catheter with a Dacron Cuff for Dialysis Short-Term Vascular Access." American Journal of Kidney Diseases, Vol. XII, Nr 6 (grudzień), 1988: strony 492-498.

Schwab, Steve J. et al. "Prospective Evaluation of a Dacron Cuffed Hemodialysis Catheter for Prolonged Use." American Journal of Kidney Diseases, Vol. XI, Nr 2 (Luty), 1988: strony 166-169

Kirkpatrick, W. G., Culpepper, R. M. & Sirmon, M.D., Frequency of Complications with Prolonged Femoral Vein Catheterization for Hemodialysis Access. Nephron: 1996; 73: strony 58-62.

Zaleski, G. X., Lorenz, J. M., Garofalo, R. S., Moscote, M. A., Rosenblum, J. D. and Leef, J. A. (1998). Experience with Tunneled Femoral Hemodialysis Catheters. American Journal of Radiology, 172: strony 493-496.

Heparin Sodium Summary of Product Characteristics

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EK OKUMA

Moss, Alvin H. et al. "Use of a Silicone Catheter with a Dacron Cuff for Dialysis Short-Term Vascular Access." American Journal of Kidney Diseases, Vol. XII, No. 6 (Aralık), 1988: sf. 492-498.

Schwab, Steve J. et al. "Prospective Evaluation of a Dacron Cuffed Hemodialysis Catheter for Prolonged Use." American Journal of Kidney Diseases, Vol. XI, No. 2 (Şubat), 1988: sf. 166-169

Kirkpatrick, W. G., Culpepper, R. M. & Sirmon, M.D., Frequency of Complications with Prolonged Femoral Vein Catheterization for Hemodialysis Access. Nephron: 1996; 73: sf. 58-62.

Zaleski, G. X., Lorenz, J. M., Garofalo, R. S., Moscote, M. A., Rosenblum, J. D. ve Leef, J. A. (1998). Experience with Tunneled Femoral Hemodialysis Catheters. American Journal of Radiology, 172: sf. 493-496.

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0123

IFU No. CR8637201-R1

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG.

™* Trademarks of its respective owner.

Other brands are trademarks of a Covidien company.

©2011 Covidien.

 Covidien Inc, 15 Hampshire Street, Mansfield, MA 02048 USA.

 Covidien Ireland Limited, IDA Business & Technology Park, Tullamore.

www.covidien.com

REV 12/2013

OPTION A

Product Brand:	Palindrome™
Product Descriptor:	H Chronic Catheter Kit Symmetrical Tip, Heparin Coating and Tal VenaTrac™ Stylet
REF #:	GNRC
GBU & COT:	Vascular Therapies/Dialysis
REV: (MM/DD/YYYY)	12/02/2013

ADDITIONAL JOB INFORMATION (if available)

Job #:	PC0215	Dims:	3.75" W x 11.00" H
Artwork/Part #:	CR8637201-R1	Metric Dims:	9.525 cm x 27.940 cm
Print Process:	Offset	Package Type:	IFU (Instructions)
Barcode Format:	N/A	File Name:	CR8637201-R1_GNRC_IFU
Dieline:	N/A	Temp ID:	P134_IFU_3.63x11_GNRC
Check here if:	<input checked="" type="checkbox"/> US Only <input type="checkbox"/> Generic Artwork	Linked to Generic Artwork: (if applicable)	

NUMBER OF COLORS

					
Black					
					

TINTS

					
---	---	---	---	---	---

TECHNICAL INFORMATION

Application/Version:	InDesign/CS5
Placed Images:	N/A
Fonts Used:	Myriad Pro

Comments:

Printer Notes:

COVIDIEN

15 Hampshire Street,
Mansfield, MA 02048

GBU-SPECIFIC INFORMATION

Graphic Designer:	Noel Gouveia

DESIGN AGENCY INFORMATION (for external use only)

Agency Name:	N/A
Agency Contact:	N/A